

Rabbit Anti-Goat IgG (H&L) FITC

Catalog No: #L35081



Package Size: #L35081-1 100 ul #L35081-2 500 ul #L35081-3 1 ml

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rabbit Anti-Goat IgG (H&L) FITC
Brief Description	Conjugated Secondary Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads.
Applications	IF, FC
Species Reactivity	G
Specificity	By immunoelectrophoresis and ELISA this antibody reacts specifically with Goat IgG. No antibody was detected against non immunoglobulin serum proteins.
Immunogen Description	Goat IgG
Conjugates	FITC
Formulation	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.05% Sodium Azide

Application Details

IF (1/250 - 1/2000), FC (1/250 - 1/2000)

Published Papers

el at., Gonadal development and molecular analysis revealed the critical window for sex differentiation, and E2 reversibility of XY-male spotted scat, *Scatophagus argus.*, , (2021)

PMID:

Note: This product is for in vitro research use only and is not intended for use in humans or animals.