

# Goat Anti-Mouse IgM mu chain - AF532

Catalog No: #L35234



Package Size: #L35234-1 100ul #L35234-2 500ul #L35234-3 1ml

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Goat Anti-Mouse IgM mu chain - AF532
Brief Description	Conjugated Secondary Antibody
Host Species	Goat
Clonality	Polyclonal
Purification	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Applications	E, IH, IF, FC
Species Reactivity	M
Specificity	By immunoelectrophoresis and ELISA this antibody reacts specifically with Mouse IgM mu chain. No antibody was detected against non immunoglobulin serum proteins.
Immunogen Description	Mouse IgM mu chain
Conjugates	AF532
Formulation	1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide

## Application Details

IH (1/100 - 1/1000), IF (1/100 - 1/1000), FC (1/1000 - 1/4000), E (Use at an assay dependent concentration)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.