## **Product Datasheet**

## Rabbit Anti-Bovine IgG (H&L) - AF568

Catalog No: #L35323

Package Size: #L35323-1 100 ul #L35323-2 500 ul #L35323-3 1 ml



Orders: order@signalwayantibody.com Support: tech@signalwayantibody.com

Description	
Product Name	Rabbit Anti-Bovine IgG (H&L) - AF568
Brief Description	Conjugated Secondary Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose
	beads.
Applications	E, IH, IF, FC
Species Reactivity	В
Specificity	By immunoelectrophoresis and ELISA this antibody reacts specifically with Bovine IgG. No antibody was
	detected against non immunoglobulin serum proteins.

1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium

## **Application Details**

Immunogen Description

Conjugates

Formulation

IH (1/100 - 1/1000), IF (1/100 - 1/1000), FC (1/1000 - 1/4000), E (Use at an assay dependent concentration)

Bovine IgG

AF568

Azide

Note: This product is for in vitro research use only and is not intended for use in humans or animals.