

# Mouse Anti-Human Ig kappa chain - AF594

Catalog No: #L35394



Package Size: #L35394-1 100 ul #L35394-2 500 ul #L35394-3 1 ml

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Anti-Human Ig kappa chain - AF594
Brief Description	Conjugated Secondary Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads.
Applications	E, IH, IF, FC
Species Reactivity	H
Specificity	By immunoelectrophoresis and ELISA this antibody reacts specifically with Human Ig kappa chain. No antibody was detected against non immunoglobulin serum proteins.
Immunogen Description	Human Ig kappa chain
Conjugates	AF594
Formulation	1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide

## Application Details

IH (1/100 - 1/1000), IF (1/100 - 1/1000), FC (1/1000 - 1/4000), E (Use at an assay dependent concentration)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.