

Mouse Anti-Human Ig kappa chain - AF633

Catalog No: #L35437



Package Size: #L35437-1 100 ul #L35437-2 500 ul #L35437-3 1 ml

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Anti-Human Ig kappa chain - AF633
Brief Description	Conjugated Secondary Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads.
Applications	E, IH, IF, FC
Species Reactivity	H
Specificity	By immunoelectrophoresis and ELISA this antibody reacts specifically with Human Ig kappa chain. No antibody was detected against non immunoglobulin serum proteins.
Immunogen Description	Human Ig kappa chain
Conjugates	AF633
Formulation	1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide

Application Details

IH (1/100 - 1/1000), IF (1/100 - 1/1000), FC (1/1000 - 1/4000), E (Use at an assay dependent concentration)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.