

# Rabbit Anti-Duck IgY (H&L) - AF647

Catalog No: #L35457



Package Size: #L35457-1 100 ul #L35457-2 500 ul #L35457-3 1 ml

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rabbit Anti-Duck IgY (H&L) - AF647
Brief Description	Conjugated Secondary Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Applications	E, IH, IF, FC
Species Reactivity	Dk
Specificity	By immunoelectrophoresis and ELISA this antibody reacts specifically with Duck IgY. No antibody was detected against non immunoglobulin serum proteins.
Immunogen Description	Duck IgY
Conjugates	AF647
Formulation	1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide

## Application Details

IH (1/100 - 1/1000), IF (1/100 - 1/1000), FC (1/1000 - 1/4000), E (Use at an assay dependent concentration)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.