

Rabbit Anti-Pig IgG (H&L) - AF647

Catalog No: #L35468



Package Size: #L35468-1 100 ul #L35468-2 500 ul #L35468-3 1 ml

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name | Rabbit Anti-Pig IgG (H&L) - AF647 |
| Brief Description | Conjugated Secondary Antibody |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads. |
| Applications | E, IH, IF, FC |
| Species Reactivity | P |
| Specificity | By immunoelectrophoresis and ELISA this antibody reacts specifically with Pig IgG. No antibody was detected against non immunoglobulin serum proteins. |
| Immunogen Description | Pig IgG |
| Conjugates | AF647 |
| Formulation | 1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide |

Application Details

IH (1/100 - 1/1000), IF (1/100 - 1/1000), FC (1/1000 - 1/4000), E (Use at an assay dependent concentration)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.