

Goat Anti-Mouse IgM mu chain - AF660

Catalog No: #L35492



Package Size: #L35492-1 100ul #L35492-2 500ul #L35492-3 1ml

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|---|
| Product Name | Goat Anti-Mouse IgM mu chain - AF660 |
| Brief Description | Conjugated Secondary Antibody |
| Host Species | Goat |
| Clonality | Polyclonal |
| Purification | The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads. |
| Applications | E, IH, IF, FC |
| Species Reactivity | M |
| Specificity | By immunoelectrophoresis and ELISA this antibody reacts specifically with Mouse IgM mu chain. No antibody was detected against non immunoglobulin serum proteins. |
| Immunogen Description | Mouse IgM mu chain |
| Conjugates | AF660 |
| Formulation | 1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide |

Application Details

IH (1/100 - 1/1000), IF (1/100 - 1/1000), FC (1/1000 - 1/4000), E (Use at an assay dependent concentration)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.