

Rabbit Anti-Guinea Pig IgG (H&L) - AF660

Catalog No: #L35502



Package Size: #L35502-1 100 ul #L35502-2 500 ul #L35502-3 1 ml

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rabbit Anti-Guinea Pig IgG (H&L) - AF660
Brief Description	Conjugated Secondary Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Applications	E, IH, IF, FC
Species Reactivity	Gp
Specificity	By immunoelectrophoresis and ELISA this antibody reacts specifically with Guinea Pig IgG. No antibody was detected against non immunoglobulin serum proteins.
Immunogen Description	Guinea Pig IgG
Conjugates	AF660
Formulation	1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide

Application Details

IH (1/100 - 1/1000), IF (1/100 - 1/1000), FC (1/1000 - 1/4000), E (Use at an assay dependent concentration)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.