

# Goat Anti-Rabbit IgG (H&L) - AF680

Catalog No: #L35537



Package Size: #L35537-1 100 ul #L35537-2 500 ul #L35537-3 1 ml

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Goat Anti-Rabbit IgG (H&L) - AF680
Brief Description	Conjugated Secondary Antibody
Host Species	Goat
Clonality	Polyclonal
Purification	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Applications	E, IH, IF, FC
Species Reactivity	Rb
Specificity	By immunoelectrophoresis and ELISA this antibody reacts specifically with Rabbit IgG. No antibody was detected against non immunoglobulin serum proteins.
Immunogen Description	Rabbit IgG
Conjugates	AF680
Formulation	1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide

## Application Details

IH (1/100 - 1/1000), IF (1/100 - 1/1000), FC (1/1000 - 1/4000), E (Use at an assay dependent concentration)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.