

# Rabbit Anti-Human IgM - AF680

Catalog No: #L35551



Package Size: #L35551-1 100 ul #L35551-2 500 ul #L35551-3 1 ml

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rabbit Anti-Human IgM - AF680
Brief Description	Conjugated Secondary Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Applications	E, IH, IF, FC
Species Reactivity	H
Specificity	By immunoelectrophoresis and ELISA this antibody reacts specifically with Human IgM. No antibody was detected against non immunoglobulin serum proteins.
Immunogen Description	Human IgM
Conjugates	AF680
Formulation	1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide

## Application Details

IH (1/100 - 1/1000), IF (1/100 - 1/1000), FC (1/1000 - 1/4000), E (Use at an assay dependent concentration)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.