

# Rabbit Anti-Pig IgG (H&L) - AF790

Catalog No: #L35640



Package Size: #L35640-1 100 ul #L35640-2 500 ul #L35640-3 1 ml

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rabbit Anti-Pig IgG (H&L) - AF790
Brief Description	Conjugated Secondary Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Applications	E, IH, IF, FC
Species Reactivity	P
Specificity	By immunoelectrophoresis and ELISA this antibody reacts specifically with Pig IgG. No antibody was detected against non immunoglobulin serum proteins.
Immunogen Description	Pig IgG
Conjugates	AF790
Formulation	1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide

## Application Details

IH (1/100 - 1/1000), IF (1/100 - 1/1000), FC (1/1000 - 1/4000), E (Use at an assay dependent concentration)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.