

Chicken Heme Oxygenase 1, Decycling (HO1) ELISA Kit

Catalog No: #EK20344

Package Size: #EK20344-1 48T #EK20344-2 96T



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Chicken Heme Oxygenase 1, Decycling (HO1) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Chicken (Gallus)
Storage	12 months

Application Details

Detect Range: 0.31-20ng/mL

Sensitivity: 0.12ng/mL

Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 uL

Assay Time: 1-4.5h

Detection wavelength: 450 nm

Product Description

Detection Method: sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.