

Chicken NADH Dehydrogenase, Quinone 1 (NQO1) ELISA Kit

Catalog No: #EK20403

Package Size: #EK20403-1 48T #EK20403-2 96T



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Chicken NADH Dehydrogenase, Quinone 1 (NQO1) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Chicken (Gallus)
Storage	12 months

Application Details

Detect Range:0.31-20ng/mL

Sensitivity: 0.14 ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 uL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.