

ZNF449 Antibody

Catalog No: #46719



Package Size: #46719-1 50ul #46719-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ZNF449 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total ZNF449 protein.
Immunogen Type	peptide
Immunogen Description	Synthetic protein corresponding to residues near the N terminal of human ZNF449
Target Name	ZNF449
Other Names	ZSCAN19
Accession No.	Swiss-Prot:Q6P9G9NCBI Gene ID:203523NCBI Protein:BC065938
Concentration	1.2mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C


Application Details

Immunohistochemistry: 1: 40-200

Images



The image on the left is immunohistochemistry of paraffin-embedded Human liver cancer tissue using 46719(ZNF449 Antibody) at dilution 1/50, on the right is treated with fusion protein. (Original magnification: x200)



The image on the left is immunohistochemistry of paraffin-embedded Human ovarian cancer tissue using 46719(ZNF449 Antibody) at dilution 1/50, on the right is treated with fusion protein. (Original magnification: x200)

Background

This gene encodes a nuclear protein that likely functions as a transcription factor. The protein includes an N-terminal SCAN domain, and seven C2H2-type zinc finger motifs.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.