

CD23 (3D1) Mouse mAb

Catalog No: #C26203

Package Size: #C26203-AF350 100ul #C26203-AF405 100ul #C26203-AF488 100ul #C26203-AF555 100ul #C26203-AF594 100ul #C26203-AF647 100ul #C26203-AF680 100ul #C26203-AF750 100ul #C26203-Biotin 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	CD23 (3D1) Mouse mAb
Host Species	Ms
Clonality	Monoclonal
Clone No.	3D1-2G4-9B6
Isotype	IgG1
Purification	Affinity purified
Species Reactivity	Hu,Ms,Rt
Immunogen Description	Synthetic Peptide of CD23
Conjugates	Unconjugated
Modification	Unmodification
Other Names	Fc fragment of IgE low affinity II receptor for, FCER 2, Lymphocyte antigen CD23, Lymphocyt
Accession No.	SwissProtID:P06734Gene ID:2208
Calculated MW	28,38,49kDa
Cell Localization	Cell membrane, Membrane, Secreted
Tag Info	The antibody detects endogenous CD23 proteins.
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	-20°C/1 year

Application Details

IF: 1:50-200 IHC: 1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.