

Ki 67 Mouse Monoclonal Conjugated Antibody

Catalog No: #C38019

Package Size: #C38019-AF350 100ul #C38019-AF405 100ul #C38019-AF488 100ul #C38019-AF555 100ul #C38019-AF594 100ul
 #C38019-AF647 100ul #C38019-AF680 100ul #C38019-AF750 100ul #C38019-Biotin 100ul #C38019-Conjugated 50ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ki 67 Mouse Monoclonal Conjugated Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	4A8
Purification	Affinity purification using immunogen.
Applications	IF
Species Reactivity	Hu
Specificity	The Ki 67 Mouse Monoclonal antibody detects endogenous Ki 67 protein.s
Target Name	Ki 67
Other Names	antigen identified by monoclonal antibody Ki-67; Antigen KI-67; KI67; KIA; MKI67
Accession No.	Swiss-Prot#:P46013
Concentration	1.0mg/ml
Formulation	Mouse IgG1 in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at 4°C in dark for 6 months

Application Details

IF:1:50-1:200

Background

KI-67 is a nuclear protein that is associated with and may be necessary for cellular proliferation. Furthermore it is associated with ribosomal RNA transcription. Inactivation of antigen KI-67 leads to inhibition of ribosomal RNA synthesis. The Ki-67 protein (also known as MKI67) is a cellular marker for proliferation. It is strictly associated with cell proliferation. During interphase, the Ki-67 antigen can be exclusively detected within the cell nucleus, whereas in mitosis most of the protein is relocated to the surface of the chromosomes. Ki-67 protein is present during all active phases of the cell cycle (G1, S, G2, and mitosis), but is absent from resting cells (G0).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.