

CD31 Polyclonal Conjugated Antibody

Catalog No: #C40699



Package Size: #C40699-AF350 100ul #C40699-AF405 100ul #C40699-AF488 100ul #C40699-AF555 100ul #C40699-AF594 100ul
 #C40699-AF647 100ul #C40699-AF680 100ul #C40699-AF750 100ul #C40699-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD31 Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Species Reactivity	Hu Ms
Specificity	CD31 Polyclonal Antibody detects endogenous levels of CD31 protein.
Immunogen Description	Synthesized peptide derived from human CD31 around the non-phosphorylation site of Y713.
Target Name	CD31
Other Names	PECAM1; Platelet endothelial cell adhesion molecule; PECAM-1; EndoCAM; GPIIA'; PECA1; CD antigen CD31
Accession No.	Swiss-Prot: P16284NCBI Gene ID: 5175
SDS-PAGE MW	82kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/20000.

Not yet tested in other applications.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.