

# c-Fos Polyclonal Conjugated Antibody

Catalog No: #C40731

Package Size: #C40731-AF350 100ul #C40731-AF405 100ul #C40731-AF488 100ul #C40731-AF555 100ul #C40731-AF594 100ul  
 #C40731-AF647 100ul #C40731-AF680 100ul #C40731-AF750 100ul #C40731-Biotin 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
 Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	c-Fos Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Species Reactivity	Hu Ms Rt
Specificity	c-Fos Polyclonal Antibody detects endogenous levels of c-Fos protein.
Immunogen Description	Synthesized peptide derived from the N-terminal region of human c-Fos.
Target Name	c-Fos
Other Names	FOS; G0S7; Proto-oncogene c-Fos; Cellular oncogene fos; G0/G1 switch regulatory protein 7
Accession No.	Swiss-Prot: P01100NCBI Gene ID: 2353
SDS-PAGE MW	62kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at 4°C in dark for 6 months

## Application Details

Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/10000.

Not yet tested in other applications.

**Note:** This product is for in vitro research use only and is not intended for use in humans or animals.