

# Nrf2 Polyclonal Conjugated Antibody

Catalog No: #C41255

Package Size: #C41255-AF350 100ul #C41255-AF405 100ul #C41255-AF488 100ul #C41255-AF555 100ul #C41255-AF594 100ul  
 #C41255-AF647 100ul #C41255-AF680 100ul #C41255-AF750 100ul #C41255-Biotin 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
 Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Nrf2 Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IF
Species Reactivity	Hu Ms Rt
Specificity	Nrf2 Polyclonal Antibody detects endogenous levels of Nrf2 protein.
Immunogen Description	Synthesized peptide derived from the C-terminal region of human Nrf2.
Target Name	Nrf2
Other Names	NFE2L2; NRF2; Nuclear factor erythroid 2-related factor 2; NF-E2-related factor 2; NFE2-related factor 2; HEBP1; Nuclear factor; erythroid derived 2, like 2
Accession No.	Swiss-Prot: Q16236NCBI Gene ID: 4780
SDS-PAGE MW	68kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/40000.

Not yet tested in other applications.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.