

PTEN Polyclonal Conjugated Antibody

Catalog No: #C41367

Package Size: #C41367-AF350 100ul #C41367-AF405 100ul #C41367-AF488 100ul #C41367-AF555 100ul #C41367-AF594 100ul
 #C41367-AF647 100ul #C41367-AF680 100ul #C41367-AF750 100ul #C41367-Biotin 100ul

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	PTEN Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Species Reactivity	Hu Ms Rt
Specificity	PTEN Polyclonal Antibody detects endogenous levels of PTEN protein.
Immunogen Description	Synthesized peptide derived from human PTEN around the non-phosphorylation site of S385.
Target Name	PTEN
Other Names	PTEN; MMAC1; TEP1; Phosphatidylinositol 3; 4,5-trisphosphate 3-phosphatase and dual-specificity protein phosphatase PTEN; Mutated in multiple advanced cancers 1; Phosphatase and tensin homolog
Accession No.	Swiss-Prot: P60484NCBI Gene ID: 5728
SDS-PAGE MW	50kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/40000.

Not yet tested in other applications.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.