

Synuclein- α Polyclonal Conjugated Antibody

Catalog No: #C41474



Package Size: #C41474-AF350 100ul #C41474-AF405 100ul #C41474-AF488 100ul #C41474-AF555 100ul #C41474-AF594 100ul
 #C41474-AF647 100ul #C41474-AF680 100ul #C41474-AF750 100ul #C41474-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Synuclein- α Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Species Reactivity	Hu Ms Rt
Specificity	Synuclein- α Polyclonal Antibody detects endogenous levels of Synuclein- α protein.
Immunogen Description	Synthesized peptide derived from human Synuclein- α around the non-phosphorylation site of Y133.
Target Name	Synuclein- α
Other Names	SNCA; NACP; PARK1; Alpha-synuclein; Non-A beta component of AD amyloid; Non-A4 component of amyloid precursor; NACP
Accession No.	Swiss-Prot: P37840NCBI Gene ID: 6622
SDS-PAGE MW	14kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/20000.

Not yet tested in other applications.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.