

# MouseB NO(NO) ELISA Kit

Catalog No: #EK20448



Package Size: #EK20448-1 48T #EK20448-2 96T

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	MouseB NO(NO) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Mouse
Storage	12 months

## Application Details

Detect Range:31.2-2000pg/mL

Sensitivity:10.1pg/mL

Sample Type:serum, plasma, tissue homogenates or other biological fluids.

Sample Volume: 100 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.