

Sheep Interferon Beta (IFN β) ELISA Kit

Catalog No: #EK20136



Package Size: #EK20136-1 48T #EK20136-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Sheep Interferon Beta (IFN β) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA Kit
Species Reactivity	Sheep
Storage	12 months

Application Details

Detect Range:31.2-2000pg/mL

Sensitivity:10.1pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 uL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.