

# Porcine Interleukin 1 Beta (IL1b) ELISA Kit

Catalog No: #EK20124



Package Size: #EK20124-1 48T #EK20124-2 96T

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Porcine Interleukin 1 Beta (IL1b) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA Kit
Storage	12 months

## Application Details

Detect Range:62.5-4000pg/mL

Sensitivity:26.2pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 uL

Assay Time:1-4.5h

Detection wavelength:450 nm

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.