

UCHL1 Rabbit Polyclonal Antibody

Catalog No: #53014

Package Size: #53014-1 50ul #53014-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	UCHL1 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human PGP9.5/PGP9.5/UCHL1 (NP_004172.2).
Other Names	HEL-117;HEL-S-53;NDGOA;PARK5;PGP 9.5;PGP9.5;PGP95;SPG79;Uch-L1;UCHL1
Accession No.	Uniprot:P09936GeneID:7345
Calculated MW	24kDa
SDS-PAGE MW	25kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using PGP9.5/PGP9.5/UCHL1 antibody.

Background

The protein encoded by this gene belongs to the peptidase C12 family. This enzyme is a thiol protease that hydrolyzes a peptide bond at the C-terminal glycine of ubiquitin. This gene is specifically expressed in the neurons and in cells of the diffuse neuroendocrine system. Mutations in this gene may be associated with Parkinson disease.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.