

NUR77 Rabbit Polyclonal Antibody

Catalog No: #53079

Package Size: #53079-1 50ul #53079-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	NUR77 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human NUR77 (NP_775180.1).
Other Names	NR4A1;GFRP1;HMR;N10;NAK-1;NGFIB;NP10;NUR77;TR3
Accession No.	Uniprot:P22736GeneID:3164
Calculated MW	33kDa/64kDa/65kDa
SDS-PAGE MW	71kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of SH-SY5Y cells, using NUR77 antibody.

Background

This gene encodes a member of the steroid-thyroid hormone-retinoid receptor superfamily. Expression is induced by phytohemagglutinin in human lymphocytes and by serum stimulation of arrested fibroblasts. The encoded protein acts as a nuclear transcription factor. Translocation of the protein from the nucleus to mitochondria induces apoptosis. Multiple transcript variants encoding different isoforms have been found for this gene.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.