

TRIM31 Rabbit Polyclonal Antibody

Catalog No: #53191



Package Size: #53191-1 50ul #53191-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	TRIM31 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human TRIM31 (NP_008959.3).
Other Names	TRIM31;C6orf13;HCG1;HCGI;RNF
Accession No.	Swiss Prot:Q9BZY9GeneID:11074
Calculated MW	30kDa/48kDa
SDS-PAGE MW	48kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB 1:500 - 1:2000IF 1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using TRIM31 at 1:1000 dilution.

Background

This gene encodes a protein that functions as an E3 ubiquitin-protein ligase. This gene shows altered expression in certain tumors and may be a negative regulator of cell growth. Alternative splicing results in multiple transcript variants.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.