

JUND Rabbit Polyclonal Antibody

Catalog No: #53290

Package Size: #53290-1 50ul #53290-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	JUND Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human JUND (NP_005345.3).
Other Names	JUND;AP-1
Accession No.	Swiss Prot:P17535GenelD:3727
Calculated MW	35kDa
SDS-PAGE MW	35kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB □ 1:500 - 1:2000 IHC □ 1:50 - 1:100

Images



Western blot analysis of extracts of various cell lines, using JUND at 1:3000 dilution.

Background

The protein encoded by this intronless gene is a member of the JUN family, and a functional component of the AP1 transcription factor complex. This protein has been proposed to protect cells from p53-dependent senescence and apoptosis. Alternative translation initiation site usage results in the production of different isoforms (PMID:12105216).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.