

NCKAP1 Rabbit Polyclonal Antibody

Catalog No: #53316

Package Size: #53316-1 50ul #53316-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	NCKAP1 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human NCKAP1 (NP_038464.1).
Other Names	NCKAP1;HEM2;NAP1;NAP125;p125Nap1
Accession No.	Swiss Prot:Q9Y2A7GeneID:10787
Calculated MW	128kDa/129kDa
SDS-PAGE MW	125kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB □ 1:500 - 1:2000 IF □ 1:50 - 1:100

Images



Western blot analysis of extracts of various cell lines, using NCKAP1 at 1:3000 dilution.

Background

-

Note: This product is for in vitro research use only and is not intended for use in humans or animals.