

VAPA Rabbit Polyclonal Antibody

Catalog No: #53470



Package Size: #53470-1 50ul #53470-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	VAPA Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human VAPA (NP_919415.2).
Other Names	VAPA;VAP-33;VAP-A;VAP33;hVAP-33
Accession No.	Swiss Prot:Q9P0L0GeneID:9218
Calculated MW	27kDa/32kDa
SDS-PAGE MW	28kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB □ 1:500 - 1:2000 IHC □ 1:50 - 1:200 IF □ 1:50 - 1:100

Images



Western blot analysis of extracts of various cell lines, using VAPA at 1:3000 dilution.



Immunohistochemistry of paraffin-embedded human stomach using VAPA at dilution of 1:100 (40x lens).

Background

The protein encoded by this gene is a type IV membrane protein. It is present in the plasma membrane and intracellular vesicles. It may also be associated with the cytoskeleton. This protein may function in vesicle trafficking, membrane fusion, protein complex assembly and cell motility. Alternative splicing occurs at this locus and two transcript variants encoding distinct isoforms have been identified.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.