

DLL4 Rabbit Polyclonal Antibody

Catalog No: #29689

Package Size: #29689-1 50ul #29689-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	DLL4 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human DLL4 (NP_061947.1).
Other Names	DLL4;AOS6;hdelta2
Accession No.	Swiss Prot:Q9NR61GeneID:54567
Calculated MW	74kDa
SDS-PAGE MW	75kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB □ 1:500 - 1:2000

IF □ 1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using DLL4 at 1:3000 dilution.

Background

This gene is a homolog of the Drosophila delta gene. The delta gene family encodes Notch ligands that are characterized by a DSL domain, EGF repeats, and a transmembrane domain.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.