

DNAJC9 Rabbit Polyclonal Antibody

Catalog No: #53520

Package Size: #53520-1 50ul #53520-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	DNAJC9 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human DNAJC9 (NP_056005.1).
Other Names	DNAJC9;HDJC9;JDD1;SB73
Accession No.	Swiss Prot:Q8WXX5GenelD:23234
Calculated MW	29kDa
SDS-PAGE MW	36kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB □ 1:500 - 1:2000 IF □ 1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using DNAJC9 at 1:3000 dilution.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.