

VPS54 Rabbit Polyclonal Antibody

Catalog No: #53527

Package Size: #53527-1 50ul #53527-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	VPS54 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human VPS54 (NP_001005739.1).
Other Names	VPS54;HCC8;PPP1R164;SLP-8p;VPS54L;WR;hVps54L
Accession No.	Swiss Prot:Q9P1Q0GeneID:51542
Calculated MW	15kDa/93kDa/109kDa/110kDa
SDS-PAGE MW	110kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB □ 1:500 - 1:2000 IF □ 1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using VPS54 at 1:3000 dilution.

Background

This gene encodes for a protein that in yeast forms part of a trimeric vacuolar-protein-sorting complex that is required for retrograde transport of proteins from prevacuoles to the late Golgi compartment. As in yeast, mammalian Vps54 proteins contain a coiled-coil region and dileucine motifs. Alternative splicing results in multiple transcript variants encoding different isoforms

Note: This product is for in vitro research use only and is not intended for use in humans or animals.