

VTN Rabbit Polyclonal Antibody

Catalog No: #53652

Package Size: #53652-1 50ul #53652-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	VTN Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human Vitronectin (NP_000629.3).
Other Names	VTN;V75;VN;VNT
Accession No.	Uniprot:P04004GeneID:7448
Calculated MW	54kDa
SDS-PAGE MW	70kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of human serum, using Vitronectin antibody.

Background

The protein encoded by this gene is a member of the pexin family. It is found in serum and tissues and promotes cell adhesion and spreading, inhibits the membrane-damaging effect of the terminal cytolytic complement pathway, and binds to several serpin serine protease inhibitors. It is a secreted protein and exists in either a single chain form or a clipped, two chain form held together by a disulfide bond.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.