

NR5A2 Rabbit Polyclonal Antibody

Catalog No: #53701

Package Size: #53701-1 50ul #53701-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	NR5A2 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human NR5A2 (NP_003813.1).
Other Names	NR5A2;B1F;B1F2;CPF;FTF;FTZ-F1;FTZ-F1beta;LRH-1;LRH1;hB1F-2
Accession No.	Swiss Prot:O00482GeneID:2494
Calculated MW	42kDa/53kDa/56kDa/61kDa
SDS-PAGE MW	61kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB □ 1:500 - 1:2000 IHC □ 1:50 - 1:200 IF □ 1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using NR5A2 at 1:1000 dilution.

Background

The protein encoded by this gene is a DNA-binding zinc finger transcription factor and is a member of the fushi tarazu factor-1 subfamily of orphan nuclear receptors. The encoded protein is involved in the expression of genes for hepatitis B virus and cholesterol biosynthesis, and may be an important regulator of embryonic development.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.