

XAB2 Rabbit Polyclonal Antibody

Catalog No: #53721

Package Size: #53721-1 50ul #53721-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	XAB2 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human XAB2 (NP_064581.2).
Other Names	XAB2;HCNP;HCRN;NTC90;SYF1
Accession No.	Swiss Prot:Q9HCS7GeneID:56949
Calculated MW	100kDa
SDS-PAGE MW	105kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB □ 1:500 - 1:2000 IF □ 1:50 - 1:100

Images



Western blot analysis of extracts of various cell lines, using XAB2 at 1:3000 dilution.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.