

DLL1 Rabbit Polyclonal Antibody

Catalog No: #53838



Package Size: #53838-1 50ul #53838-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	DLL1 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human DLL1 (NP_005609.3).
Other Names	DLL1;DELTA1;DL1;Delta
Accession No.	Swiss Prot:O00548GeneID:28514
Calculated MW	27kDa/78kDa
SDS-PAGE MW	78kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB □ 1:500 - 1:2000 IHC □ 1:50 - 1:200


Images



Western blot analysis of extracts of various cell lines, using DLL1 at 1:3000 dilution.



Immunohistochemistry of paraffin-embedded rat pancreas using DLL1 at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded human breast cancer using DLL1 at dilution of 1:100 (40x lens).

Background

DLL1 is a human homolog of the Notch Delta ligand and is a member of the delta/serrate/jagged family. It plays a role in mediating cell fate decisions during hematopoiesis. It may play a role in cell-to-cell communication.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.