

## DDX59 Rabbit Polyclonal Antibody

Catalog No: #53948

Package Size: #53948-1 50ul #53948-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	DDX59 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Rat
Immunogen Description	Recombinant fusion protein of human DDX59 (NP_001026895.2).
Other Names	DDX59;OFD5;ZNHIT5
Accession No.	Swiss Prot:Q5T1V6GeneID:83479
Calculated MW	64kDa/68kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

## Application Details

WB □ 1:500 - 1:2000 IHC □ 1:50 - 1:200

## Images



Immunohistochemistry of paraffin-embedded rat testis using DDX59 at dilution of 1:100 (40x lens).

## Background

-

Note: This product is for in vitro research use only and is not intended for use in humans or animals.