

LDHAL6B Rabbit Polyclonal Antibody

Catalog No: #53950

Package Size: #53950-1 50ul #53950-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	LDHAL6B Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Mouse,Rat
Immunogen Description	Recombinant fusion protein of human LDHAL6B (NP_149972.1).
Other Names	LDHAL6B;LDH6B;LDHAL6;LDHL
Accession No.	Swiss Prot:Q9BYZ2GeneID:92483
Calculated MW	41kDa
SDS-PAGE MW	42kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB □ 1:500 - 1:2000 IF □ 1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using LDHAL6B at 1:1000 dilution.



Immunofluorescence analysis of rat testis using LDHAL6B at dilution of 1:100. Blue: DAPI for nuclear staining.

Background

-

Note: This product is for in vitro research use only and is not intended for use in humans or animals.