

OCIAD2 Rabbit Polyclonal Antibody

Catalog No: #53953

Package Size: #53953-1 50ul #53953-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	OCIAD2 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human OCIAD2 (NP_689611.1).
Other Names	OCIAD2
Accession No.	Swiss Prot:Q56VL3GeneID:132299
Calculated MW	11kDa/16kDa
SDS-PAGE MW	17kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB □ 1:500 - 1:2000 IF □ 1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using OCIAD2 at 1:1000 dilution.

Background

-

Note: This product is for in vitro research use only and is not intended for use in humans or animals.