

## VILL Rabbit Polyclonal Antibody

Catalog No: #54037

Package Size: #54037-1 50ul #54037-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	VILL Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human VILL (NP_056957.3).
Other Names	VILL
Accession No.	Swiss Prot:O15195GeneID:50853
Calculated MW	94kDa/95kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

## Application Details

WB □ 1:200 - 1:2000 IF □ 1:50 - 1:200

## Images



Immunofluorescence analysis of L929 cells using VILL at dilution of 1:100. Blue: DAPI for nuclear staining.

## Background

The protein encoded by this gene belongs to the villin/gelsolin family. It contains 6 gelsolin-like repeats and a headpiece domain. It may play a role in actin-bundling.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.