

DNER Rabbit Polyclonal Antibody

Catalog No: #54052



Package Size: #54052-1 50ul #54052-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	DNER Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human DNER (NP_620711.3).
Other Names	DNER;UNQ26;bet
Accession No.	Swiss Prot:Q8NFT8GeneID:92737
Calculated MW	78kDa
SDS-PAGE MW	80kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB □ 1:200 - 1:2000 IF □ 1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using DNER at 1:1000 dilution.



Immunofluorescence analysis of L929 cells using DNER at dilution of 1:100. Blue: DAPI for nuclear staining.

Background

-

Note: This product is for in vitro research use only and is not intended for use in humans or animals.