

ZNF263 Rabbit Polyclonal Antibody

Catalog No: #54227



Package Size: #54227-1 50ul #54227-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ZNF263 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human ZNF263 (NP_005732.2).
Other Names	ZNF263;FPM315;ZKSCAN12;ZSCAN44
Accession No.	Swiss Prot:O14978GeneID:10127
Calculated MW	77kDa
SDS-PAGE MW	77kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB □ 1:2000 - 1:5000 IHC □ 1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using ZNF263 at 1:1000 dilution.

Background

-

Note: This product is for in vitro research use only and is not intended for use in humans or animals.