

NGEF Rabbit Polyclonal Antibody

Catalog No: #54231

Package Size: #54231-1 50ul #54231-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	NGEF Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human NGEF (NP_062824.2).
Other Names	NGEF;ARHGEF27;EPHEXIN;ephexin-1
Accession No.	Swiss Prot:Q8N5V2GeneID:25791
Calculated MW	26kDa/71kDa/82kDa
SDS-PAGE MW	82kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB □ 1:500 - 1:2000 IHC □ 1:50 - 1:100 IF □ 1:50 - 1:100

Images



Western blot analysis of extracts of various cell lines, using NGEF at 1:1000 dilution.

Background

-

Note: This product is for in vitro research use only and is not intended for use in humans or animals.