

OTUD7B Rabbit Polyclonal Antibody

Catalog No: #54409



Package Size: #54409-1 50ul #54409-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	OTUD7B Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human OTUD7B (NP_064590.2).
Other Names	OTUD7B;CEZANNE;ZA20D1
Accession No.	Swiss Prot:Q6GQQ9GeneID:56957
Calculated MW	88kDa/92kDa
SDS-PAGE MW	95Da
Formulation	Buffer: PBS with 0.02% sodium azide, pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB □ 1:500 - 1:2000 IHC □ 1:50 - 1:200 IF □ 1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using OTUD7B at 1:1000 dilution.



Immunofluorescence analysis of L929 cells using OTUD7B at dilution of 1:100. Blue: DAPI for nuclear staining.

Background

-

Note: This product is for in vitro research use only and is not intended for use in humans or animals.