

NEK9 Rabbit Polyclonal Antibody

Catalog No: #54465

Package Size: #54465-1 50ul #54465-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	NEK9 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human NEK9 (NP_149107.4).
Other Names	APUG;LCCS10;NC;NERCC;NERCC1;NEK9
Accession No.	Swiss Prot:Q8TD19GeneID:91754
Calculated MW	107kDa
SDS-PAGE MW	107kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB □ 1:2000 - 1:5000 IHC □ 1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using NEK9 at 1:1000 dilution.

Background

This gene encodes a member of the NimA (never in mitosis A) family of serine/threonine protein kinases. The encoded protein is activated in mitosis and, in turn, activates other family members during mitosis. This protein also mediates cellular processes that are essential for interphase progression. [provided by RefSeq, Jul 2016]

Note: This product is for in vitro research use only and is not intended for use in humans or animals.