

## MT-ATP6 Rabbit Polyclonal Antibody

Catalog No: #54482



Package Size: #54482-1 50ul #54482-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	MT-ATP6 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of human MT-ATP6.
Other Names	ATPase6;MTATP6;MT-ATP6
Accession No.	Swiss Prot:P00846GenelD:4508
Calculated MW	25kDa
SDS-PAGE MW	25kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


## Application Details

WB □ 1:500 - 1:2000

IHC □ 1:50 - 1:100


IF □ 1:50 - 1:100

## Images



Immunofluorescence analysis of L929 cells using MT-ATP6 Polyclonal at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.

Western blot analysis of extracts of various cell lines, using MT-ATP6 at 1:1000 dilution.



## Background

-

Note: This product is for in vitro research use only and is not intended for use in humans or animals.