

DGCR8 Rabbit Polyclonal Antibody

Catalog No: #54579

Package Size: #54579-1 50ul #54579-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	DGCR8 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of human DGCR8.
Other Names	—
Accession No.	Swiss Prot:Q8WYQ5GeneID:54487
Calculated MW	—
SDS-PAGE MW	110kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB □ 1:2000 - 1:5000 IHC □ 1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using DGCR8 at 1:1000 dilution.

Background

—

Note: This product is for in vitro research use only and is not intended for use in humans or animals.