

TACC3 Rabbit Polyclonal Antibody

Catalog No: #54582

Package Size: #54582-1 50ul #54582-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TACC3 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human
Immunogen Description	Recombinant protein of human TACC3.
Other Names	–
Accession No.	Swiss Prot:Q9Y6A5GeneID:10460
Calculated MW	–
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB □ 1:500 - 1:2000 IF □ 1:50 - 1:200

Images



Immunofluorescence analysis of U2OS cells using TACC3 at dilution of 1:100. Blue: DAPI for nuclear staining.

Background

–

Note: This product is for in vitro research use only and is not intended for use in humans or animals.